



MacGregor<sup>®</sup>  
Healthcare Ltd

# The **Qufora**<sup>®</sup> IrriSedo Rectal/Transanal irrigation range



**Qufora**<sup>®</sup> IrriSedo Mini



**Qufora**<sup>®</sup> IrriSedo Cone



**Qufora**<sup>®</sup> IrriSedo Klick



**Qufora**<sup>®</sup> IrriSedo Bed

## Product Choice Guide

**Distributed by**  
MacGregor Healthcare Ltd  
[www.macgregorhealthcare.com](http://www.macgregorhealthcare.com)  
[info@macgregorhealthcare.com](mailto:info@macgregorhealthcare.com)  
Tel: 0845 519 2040

MBH | **qufora**<sup>®</sup>





## Qufora® IrriSedo Klick



Qufora® IrriSedo Klick is a higher volume rectal irrigation with a water inflated balloon catheter to create a water seal. The system is designed for easy assembly and usage with a click connection.

The catheter has a hydrophilic coating, by dipping the catheter in water, the coating is activated and ready to use. The control unit follows a simple 1,2,3,4 sequence.

For assessment, teaching, regime and follow up please refer to the **'Decision guide'**.

- ✓ Higher volume irrigation
- ✓ Good/poor dexterity
- ✓ Can sit on or over toilet with support
- ✓ Catheter self-supporting when balloon inflated
- ✓ Choice of gravity or pump

Urge faecal incontinence/urgency

Evacuation difficulties

Constipation (slow transit/idiopathic/opioid induced/IBS-C

Neurogenic (spinal cord injury/upper motor neurone/ MS/parkinson's/spina bifida/cauda equina



## Qufora® IrriSedo Bed



The Qufora® IrriSedo Bed system is a unique way of performing rectal irrigation while the patient remains in bed.

The system instills water into the rectum via a catheter which is attached to a collection bag. This allows for controlled, hygienic irrigation with minimal odour.

For assessment, teaching, regime and follow up please refer to the **'Decision guide'**.

- ✓ Higher volume irrigation
- ✓ Rectal irrigation with patient in bed
- ✓ Carer performs irrigation
- ✓ Choice of gravity or pump

Bed bowel management

Poor balance/unable to transfer

Spinal cord injury with upper motor neurone lesion/trunk balance e.g. upper motor neurone lesion

**The above conditions are intended as a guide only.**

All patients should be individually assessed for their suitability to use transanal irrigation.

# Ordering



## Qufora® IrriSedo Mini

DESCRIPTION	CONTENTS	ORDER CODE	SUGGESTED MONTHLY ORDERING
Mini system (30 Irrigations)	1 Hand pump, 1 Cone bag, 30 Hydrophilic cones	53601-030	<b>When irrigating daily:</b> 1 Mini system (30 irrigations) if extension tubes being used 2 sets of extension tubes  <b>When irrigating every other day:</b> 1 Mini system (15 irrigations) if extension tubes being used 1 sets of extension tubes
Mini system (15 Irrigations)	1 Hand pump, 1 Cone bag, 15 Hydrophilic cones	53601-015	
Extension tube (15 Irrigation set)	15 Extension tubes	52715-015	



## Qufora® IrriSedo Cone

DESCRIPTION	CONTENTS	ORDER CODE	SUGGESTED MONTHLY ORDERING
Cone starter set (5 Irrigations)	1 Water bag with pump, 1 Cone bag, 1 Hook, 5 Hydrophilic cones	QTS	<b>When irrigating daily:</b> 1 cone system (15 irrigations) 1 accessory set  <b>When irrigating every other day:</b> 1 cone system (15 irrigations)
Cone system (15 Irrigations)	1 Water bag with pump, 1 Cone bag, 15 Hydrophilic cones	QTM	
Cone accessory set (15 irrigations)	15 Hydrophilic cones	52601-015	



## Qufora® IrriSedo Klick

DESCRIPTION	CONTENTS	ORDER CODE	SUGGESTED MONTHLY ORDERING
Klick starter set (5 Irrigations)	5 Catheters, 1 Klick tube, 1 Control unit with pump and velcro strap, 1 Water bag, 1 Hook, 1 Toilet bag	QKSTART	<b>First month order</b> 1 starter set 2 catheter set  <b>When irrigating daily:</b> 2 catheter sets per month 1 control unit every 3 months  <b>When irrigating every other day:</b> 1 catheter set per month 1 control unit every 6 months
Klick catheter set (15 Irrigations)	15 Catheters, 1 Klick tube, 1 Water bag	QKCATH	
Klick control unit (90 Irrigations)	1 Control unit with pump and velcro strap	QKCONTROL	



## Qufora® IrriSedo Bed

DESCRIPTION	CONTENTS	ORDER CODE	SUGGESTED MONTHLY ORDERING	
5 irrigation set	Bed system (5 Irrigations) – Ø12mm	1 Water bag with pump, 5 Rectal catheters Ø12mm, 5 Gel sachets, 5 Collection bags	53201-055	<b>When irrigating daily:</b> 2 bed systems (15 irrigations)  <b>When irrigating every other day:</b> 1 bed systems (15 irrigations)
	Bed system (5 Irrigations) – Ø21mm	1 Water bag with pump, 5 Rectal catheters Ø21mm, 5 Gel sachets, 5 Collection bags	53001-055	
	Bed system (5 Irrigations) – Ø30mm	1 Water bag with pump, 5 Rectal catheters Ø30mm, 5 Gel sachets, 5 Collection bags	53101-055	
15 irrigation set	Bed system (15 Irrigations) – Ø12mm	1 Water bag with pump, 15 Rectal catheters Ø12mm, 5 Gel sachets, 15 Collection bags	53201-015	
	Bed system (15 Irrigations) – Ø21mm	1 Water bag with pump, 15 Rectal catheters Ø21mm, 5 Gel sachets, 15 Collection bags	53001-015	
	Bed system (15 Irrigations) – Ø30mm	1 Water bag with pump, 15 Rectal catheters Ø30mm, 5 Gel sachets, 15 Collection bags	53101-015	

200114

